

EP-TRACER® 2

**schwarzer
cardiotek**
SMART. PRECISE. RELIABLE.


The EP recording system – Measure the Difference

CSP Module
available

Reliable results for your electrophysiology investigations

- Excellent signal quality
- High quality integrated stimulator with two outputs
- Connectivity to Cryo, PFA and RF generator
- 3 different view modes: real time, trigger mode, view mode
- Easy handling with one touch instructions
- Invasive blood pressure channels
- 20, 52 or 84 bipolar intracardiac channels
- CSP module with LVAT measurement, Current of Injury visualization
- 12 lead ECG
- Export of calibrated signals (in mV)
- Review possibility
- Screenshots, recordings, log files and reports
- DICOM interface

Available in **three different** versions



Smart Keyboard

Control at the palm of your hand



Service & Support

Benefit from optional service contracts for **maintenance, updates, and fast technical support** - tailored to your needs.

Contact us

Email info@schwarzercardiotek.com
Phone +49 (0)7131 2774 500